


Elefant[®]

Disposable suction/irrigation device
for laparoscopic surgery



References:

1. Abu-Rafea B., Vilos GA., Al-Obeed O., AlSheikh A., Vilos AG., MD, Al-Mandeeel M., Monopolar Electrosurgery through Single-Port Laparoscopy: A Potential Hidden Hazard for Bowel Burns, *J Minim Invasive Gynecol.* Nov-Dec 2011;18(6):
2. Alkatout I., Schollmeyer T. Nusrat A. Hawaldar, MS., Sharma N., Mettler, L., Principles and Safety Measures of Electrosurgery in Laparoscopy, *JLS* (2012)16:130-139
3. Borin J. (2011) M.C. Ost (ed.), *Robotic and Laparoscopic Reconstructive Surgery in Children and Adults*, 3 *Current Clinical Urology*, DOI 10.1007/978-1-60327-914-7_1

Customy Care / Continance Care / Wound & Skin Care / Interventional Urology

Elefant®

Disposable suction/irrigation device for laparoscopic surgery

Yankauer tip with 6 distal holes

Prevents suction of the neighbouring tissues, preserves the pneumoperitoneum.



"Non-conductive" fibreglass cannula

- Non-traumatic
- Designed to address the complications of electrothermal injuries^{1,2}
- 5 mm or 10 mm diameters
- 35 cm or 45 cm lengths

Sterile & ready to use

Single use cannula avoids cross-contamination

Ergonomic design handle

- Single-hand control
- For right or left-handed users
- Lightweight cannula: around 100 gr



2 separate irrigation lines

To connect 2 irrigation bags (Ref. no. ASP190 & ASP195)



Tubing

- Color-coded
- PVC

Reference overview

for Elefant® disposable suction/irrigation device




References for robotic surgery³

Elefant® single-use suction/irrigation device Without tubing

Ref. no.	Cannula length / cm	Cannula Ø / mm	Box of
ASP145	35	5	10
ASP150	35	10	10
ASP180	45	5	10

Elefant® single-use suction/irrigation device With double PVC tubing

Ref. no.	Cannula length / cm	Cannula Ø / mm	Tubing length / m	Tubing external Ø / mm	Tubing inner Ø / mm	Number of Irrigation lines	Box of
ASP165	35	5	3	9,5	6,4	1	5
ASP170	35	10	3	9,5	6,4	1	5
ASP185	45	5	3	9,5	6,4	1	5
ASP190	35	5	3	9,5	6,4	2	5
ASP195	45	5	3	9,5	6,4	2	5



See an Elefant® in action at www.elefantintheroom.com

